
# FISA TEHNICA PRODUS

**Denumire Produs:** Test Rapid Boson Professional  
**Producator:** Xiamen Boson Biotech Co., Ltd.

## COVID-19

Rapid SARS-CoV-2 Antigen Test Card

Using a nasopharyngeal, oropharyngeal or nasal swab to get sample fluid, antigen tests can produce results in minutes. Because these tests are faster and less expensive than molecular tests are, experts consider antigen tests more practical to use for large numbers of people. Since the antigen test can be high specificity, a positive antigen test result is considered very accurate.

- **Fast, Easy and Convenient**

Results in 15-20 minutes, simple to use  
Requires no instruments and special skills


- **Storage**

Room temperature storage & 18-month shelf life


- **Swabs Test**

Nasopharyngeal/Oropharyngeal/Nasal swabs can be used

### Test Procedure

### Interpretation of Results

**BOSON**

## Declaration of Conformity

**Manufacturer** Xiamen Boson Biotech CO., Ltd.  
90-94 Tianfeng Road, Jimei North Industrial Park,  
Xiamen, Fujian 361021, P. R. China.

**European Representative** Lotus NL B.V.  
Koningin Julianaplein 10, 1e Verd, 2595AA,  
The Hague, Netherlands.

**Device/s** Rapid SARS-CoV-2 Antigen Test Card

**Classification** Others

**Confirmative Assessment Route** 98/79/EC IVDD Annex III

We, Xiamen Boson Biotech Co., Ltd. declare that the above mentioned devices conform to the relevant provisions of the EC Council Directive 98/79 and is in accordance with the Annex III, ISO 13485:2016 Quality Management System, as implemented by the European Union's Medical Devices Regulations and the Federal and Local Authorities.

**Place, Date of Issue** Xiamen, 2020-08-04

**Signature** 

(Signed By Boson Representative)

Name: Changgong Zhang  
Title: General Manager

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Pagina 1 van 2

## In Vitro Diagnostic Reagent

EUROPEAN COMMISSION  
DIRECTORATE-GENERAL FOR HEALTH AND FOOD SAFETY  
Public health, country knowledge, crisis management  
**Health Security and Vaccination**

# Clinical Evaluation Report

**Product Name:** Rapid SARS-CoV-2 Antigen Test Card  
**Manufacturer:** Xiamen Boson Biotech Co., Ltd.

**Clinical Trial Institution:** Chongqing Public Health Medical Center  
Ankara Dr. Nafiz Korez Sincan State Hospital  
The Catholic University of Korea Uijeongbu ST. Mary's Hospital

**Person Responsible for Clinical Statistics:** Kesai Liu

**Applicant Institution:** Xiamen Boson Biotech Co., Ltd.  
**Address:** 90-94 Tianfeng Road, Jimei North Industrial Park, Xiamen  
**Contact:** Zhijuan Jia

**Report Date:** March 10, 2021

### EU health preparedness:

**A common list of COVID-19 rapid antigen tests, including those whose test results are mutually recognised, and a common standardised set of data to be included in COVID-19 test result certificates**

Agreed by the Health Security Committee  
on 17 February 2021

Xiamen Boson Biotech Co	Rapid SARS-CoV-2 Antigen Test card	Yes	Not specified		93.8% sensitivity 100% specificity NP swab			BE, BG, DE <sup>[2]</sup> , FR	CH	<a href="#">DE</a>		Yes	Yes
-------------------------	------------------------------------	-----	---------------	--	--	--	--	-----------------------------------	----	--------------------	--	-----	-----

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**COVID-19 Medical Device  
Authorization with Conditions for  
Importation or Sale Pursuant to  
Section 7 of the Interim Order**

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Pursuant to section 5 of the Interim Order Respecting the Importation and Sale of Medical Devices for Use in Relation to COVID-19, made by the Minister of Health on March 18, 2020, the medical device listed below is now authorized for sale or importation in Canada.

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Each shipment of a COVID-19 medical device that is imported into Canada must be accompanied by a copy of this authorization document. **Please ensure to highlight the Authorization reference number during the import declaration process to facilitate port entry without any delays.**

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This authorization is only valid for so long as the Interim Order Respecting the Importation and Sale of Medical Devices for Use in Relation to COVID-19 is in effect.

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XIAMEN BOSON BIOTECH CO., LTD.