FISA TEHNICA PRODUS

Denumire Produs: Test Rapid Antigen SARS-CoV-2 Realy Tech Lateral Nazal 1 Buc
Producator: Hangzhou Realy Tech Co., Ltd.

1. DESCRIERE

Test Rapid Antigen COVID Lateral Nazal (1 buc)

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Manufacturer:

Name: Hangzhou Realy Tech Co., Ltd.

Address: #2 Building, No. 763, Yuansha Village, Xinjie Street, Xiaoshan District,
311200 Hangzhou City, Zhejiang Province, PEOPLE'S REPUBLIC OF CHINA

Product/s	Catalogue number
Novel Coronavirus (SARS-Cov-2) Antigen Rapid Test Device (nasal swab)	K601416D

Category: Other Devices (All devices except Annex II and self-testing devices)

Conformity assessment route: Annex III (except Point 6) of the Directive

Applicable Standards: EN ISO 13485:2016; EN ISO 15223-1:2016; EN ISO 13612:2002/AC:2002, EN ISO 17511:2003, EN ISO 18113-1:2011, EN ISO 18113-2:2011, EN ISO 23640:2015; EN ISO 13641:2002; EN ISO 14971:2019, EN 62366-1:2015.

We, the Manufacturer, herewith declare with sole responsibility that our product/s mentioned above meet/s the provisions of the Directive 98/79/EC of the European Parliament and of the Council on In-Vitro Diagnostic Medical Devices.

We hereby explicitly appoint Luxus Lebenswelt GmbH, located at Kochstr.1,47877, Willich, Germany to act as our European Authorised Representative as defined in the aforementioned Directive.

Hangzhou 2021.3.2

(Place and date of issue)

(Signature and position)

Signed for and on behalf of the manufacturer

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Formblatt für In-vitro-Diagnostika / Form for In Vitro Diagnostic Medical Devices

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Clinical Evaluation Report

Introduction

SARS-CoV-2 is the causative agent of coronavirus disease 2019 (COVID-19), whose principal clinical manifestation is a bilateral pneumonia. The novel coronavirus belongs to the *Coronaviridae* family, genus *Betacoronavirus*, subgenus *Sarbecovirus*. Due to the efficient human-to-human transmission, the virus spread rapidly all over the world and on 11 March 2020, WHO declared COVID-19 as a pandemic. To date, SARS-CoV-2 infected more than one-hundred million people and caused the death of more than 3 million people. The virus can be isolated from both the upper and lower respiratory tract. Real-time PCR is considered the gold standard for the detection of SARS-CoV-2; however, it requires skilled personnel, is time consuming and suffer of reagents shortage. On the other hand, antigen testing is quick and can identify infectious individuals, which are those that carry a high viral load in their upper respiratory tract. Although, less sensitive than molecular test, antigen test is useful for identifying infected individuals in the first days of infection when the viral load in the upper respiratory tract is highest. Here, we tested the sensitivity and specificity of the Novel Coronavirus (SARS-CoV-2) Antigen Rapid Test using nasal swabs already tested by real-time PCR with the Allplex™-SARS-CoV-2 assay (Seegene, Seoul, South Korea).

Material and methods

Sample size

Two hundred nasal swabs collected at the COVID center of the University Hospital Tor Vergata and tested for the presence of SARS-CoV-2 by real-time PCR were included in the study. The selected samples included 100 samples positive for SARS-CoV-2 and 100 samples negative to SARS-CoV-2. All positive selected samples were amplified at a Ct value ≤ 25 . Samples were collected on May 2021.

SARS-CoV-2 real-time PCR

Nasal swabs were tested by Allplex™-SARS-CoV-2 assay (Seegene, Seoul, South Korea). The assay detects the E gene, common to the Sarbecovirus subgenus, and three specific SARS-CoV-2 genes: N, RdRp and S.

Antigen test

Samples selected on the basis of real-time PCR results were subsequently tested by Novel Coronavirus (SARS-CoV-2) Antigen Rapid Test cassette (Hangzhou Realy Tech Co., Ltd., China). Eighty μ l of nasal swab dispersed in universal transport medium (UTM) were mixed with 80 μ l of extraction buffer of the Antigen assay. Then, 80 μ l were loaded on the cassette according to the manufacturer's instruction. The appearance of control and antigen lines was checked after 1520 min as indicated by the manufacturer.

Results

Of the one-hundred positive samples with a Ct value ≤ 25 , 96 tested positive also with the Antigen assay, while 4 were negative, Table 1. The four negative samples had a mean Ct value of 22.5. The 100 negative samples by real-time PCR were all confirmed negative by the Antigen assay with a specificity of 100%, Table 2.

Conclusions

The Novel Coronavirus (SARS-Cov-2) Antigen Rapid Test showed 96% sensitivity (considering only RTPCR positive sample with a cycle threshold ≤ 25) and 100% specificity. Based on these results, it can be stated that the Novel Coronavirus (SARS-Cov-2) Antigen Rapid Test is suitable for the diagnosis of SARS-CoV-2 infection with a Detection Limit corresponding to or less than 25 CT at the molecular test (RTPCR).

Table 1. Demographics of 100 SARS-CoV-2 positive patients detected by real-time PCR (Allplex™ SARS-CoV-2 assay, Seegene) carried out on nasal swabs collected on May 2021.

NO.	Gender	Age	Experimental reagent Assessment test results		PCR test results		
			Sample type	Determination	Sample type	Determination	Ct value
1	M	18	nasal swab	Positive	nasal swab	Positive	20
2	F	57	nasal swab	Negative	nasal swab	Positive	22
3	M	18	nasal swab	Positive	nasal swab	Positive	21
4	M	20	nasal swab	Positive	nasal swab	Positive	16
5	F	18	nasal swab	Positive	nasal swab	Positive	13
6	M	26	nasal swab	Positive	nasal swab	Positive	21
7	M	80	nasal swab	Positive	nasal swab	Positive	20
8	M	37	nasal swab	Positive	nasal swab	Positive	19
9	M	26	nasal swab	Positive	nasal swab	Positive	18
10	M	17	nasal swab	Positive	nasal swab	Positive	16
11	M	41	nasal swab	Positive	nasal swab	Positive	22
12	M	45	nasal swab	Positive	nasal swab	Positive	19
13	F	46	nasal swab	Positive	nasal swab	Positive	21
14	M	38	nasal swab	Positive	nasal swab	Positive	21
15	F	37	nasal swab	Positive	nasal swab	Positive	25
16	M	37	nasal swab	Positive	nasal swab	Positive	23
17	M	68	nasal swab	Positive	nasal swab	Positive	21
18	M	43	nasal swab	Positive	nasal swab	Positive	23
19	M	56	nasal swab	Negative	nasal swab	Positive	23
20	F	48	nasal swab	Positive	nasal swab	Positive	18
21	M	19	nasal swab	Positive	nasal swab	Positive	21
22	M	46	nasal swab	Positive	nasal swab	Positive	20
23	F	51	nasal swab	Positive	nasal swab	Positive	18
24	M	72	nasal swab	Positive	nasal swab	Positive	18
25	M	14	nasal swab	Positive	nasal swab	Positive	17
26	M	34	nasal swab	Positive	nasal swab	Positive	22
27	M	61	nasal swab	Positive	nasal swab	Positive	22
28	F	33	nasal swab	Positive	nasal swab	Positive	24
29	M	19	nasal swab	Positive	nasal swab	Positive	19
30	M	51	nasal swab	Positive	nasal swab	Positive	18
31	F	46	nasal swab	Positive	nasal swab	Positive	19
32	F	78	nasal swab	Positive	nasal swab	Positive	24

Va multumim!

comenzi@medplaza.ro
office@medplaza.ro

www.medplaza.ro
0316.30.66.90